



NIRx Medical Technologies, LLC

Enhancing  
New Dimensions in Neuroimaging



**NIRSport**

[nirx.net](http://nirx.net)

## NIRSport System Description

The NIRSport is a light-weight, freely configurable, multi-channel fNIRS imaging system that combines LED illumination with active detection technology for a truly wearable brain imaging solution. These novel product features allow for a wide field of innovative applications where portability, wearability, and a small footprint are essential. This system allows for noninvasive realtime hemoglobin measurements of the cerebral cortex.



## Wearable System Solution

"The NIRSport offers 8 sources and 8 detectors (16 sources / 16 detectors in tandem mode) with a diverse array of available headgear and optical probes. Turn-key solutions are available for standard topography, neonatal/infant imaging, as well as multimodal acquisition with EEG, MRI, TMS, tDCS, and Eyetracking. NIRx offers both standardized probe positioning and custom-configurable setups.


This system has several digital input/output options for precise event marker triggering. A real-time data streaming option is available for BCI/neurofeedback applications. The NIRSport employs an open data format that is compatible with a variety of open-source analysis solutions as well as NIRx's proprietary NAVI software."




## Applications

- |                                       |                            |
|---------------------------------------|----------------------------|
| Autism                                | Motor Tasks                |
| Brain Computer Interface              | Neonatal-Infant Monitoring |
| Hyperscanning (Multi-subject Studies) | Psychiatric Disorders      |
| Infant Monitoring                     | Speech Recognition         |
| Intra Operative Monitoring            | Stroke and Rehabilitation  |
| Language                              | Traumatic Brain Injury     |
| Learning and Attention                | Wearable Studies           |

## Probe Setup



## Data Event Finder



Examples for various Advanced Plot styles during Scan

### **NIRStar Software Description**

The NIRStar software package provides a user friendly Graphical User Interface for system control (calibration and probe setup), patient monitoring, real-time cortical 2d and 3d display capabilities, and a module for hyper-scanning.

### **NIRStar Instrument Software features:**

Scalable sensor configuration

Real time event recording-display

Programmable illumination configuration

Online system performance checks

Time averaged event update display

Anatomically configurable optode arrangements: graphical 3-d digitizer

Real-time 2d, scalp, and cortical mapping with hyper-scanning capability of oxygenated, deoxygenated, and total hemoglobin



Hardware		Select System Features
NIRSport Dimensions	105 mm x 170 mm x 40 mm	2d topography • Measurement type: Change in oxy-hemoglobin, deoxy-hemoglobin and total hemoglobin concentrations • Light-weight, wearable system • Powered via laptop or tablet PC • LED Illumination: Greater temporal stability than laser source • Digital Detection, Digital Trigger I/O • Flexible imaging arrays supports multi-distance measurements • Multi-modal capability with: TMS, tDCS, MRI, EEG, Eyetracking • Real-time display, Real time data-streaming • Equipped with NAVI 2d-3d data analysis software from NIRx • Data exportable in ASCII format: compatible with other analysis software • Tandem capability: 2 separate systems synchronized to operate simultaneously • Highest dynamic range commercially available
Net weight	660 g	
Illumination type	LED	
Number of Illumination Sources	8 (16 in tandem mode)	
Detection Type	Active detection sensor (no optical fibers)	
Number of Detectors	8 (16 in tandem mode)	
Number of Channels	Up to 64 (up to 128 tandem mode)	
System scan rate	2.5Hz to 62.5Hz	
Dual Wavelength	760nm, 850nm (available up to 4 wavelengths upon request)	
Dynamic range	60 dBopt	
Sensor type	Si Photodiode	
Sensitivity	< 1 pW	
Host connection	USB 2.0 or 3.0 + USB 2.0 or 3.0 power via Host PC	
Experiment timing	TTL/CMOS (8 In/ 8 Out max.)	
E-Prime compatible	255 Input-Output Conditions	
Supply voltage	15 - 21 VDC	
Power consumption	39 VA max	
Phase type	Single Phase	
Mode of operation	Continuous Wave	
Temperature range	10 to 40 degrees C (operating), -15 to 70 degrees C (storage and transport).	
Humidity	Humidity 20 - 80% relative humidity non-condensing.	
Standard measuring cap	128 optode positions	
Static phantom	For system calibration	
Standard accessories	System carrying case, system base-plate, NIRSport backpack for wearable studies, Adapter cable to connect parallel-cable (25-pin) to 10-pin IDC trigger input	
PC Requirements	Dell Latitude E6430 with 3rd gen core i3 processor, 4 GB, 320 GB hdd, USB 2, Win 7 32-bit or equivalent configuration	

### Options

- Software Development Kit
- Custom configurable optical fibers
- Custom LED-Optical fiber coupling solution
- Integrated NIRS-EEG Head Gear
- Custom configured headgear arrays
- For real time data spooling
- Available up to 10 meters
- For multimodal studies: fNIRS-TMS, fNIRS-MRI
- With most commercial EEG systems
- Patches, Prefrontal, Whole head