

# Luna EMG



## Luna EMG: robotic device for upper and lower limb neurological and early stage orthopedic rehabilitation.

Improve a patient's muscle strength, range of motion and coordination, after strokes, SCIs, MS, TBI, orthopaedic injuries surgeries (knee/hip replacements, ankle alignment).

### Multi-Functional Isokinetic System for Training and Diagnostics



The most unique feature - **EMG-assisted movement.\*** - active training, even for very weak patients (MMT 1)

\*EMG = electromyography, electrical activity of the muscle. Robot detects the patient's muscle activity via sensors and converts it into intended movement.



This is a medical device. For your safety, use it in accordance with the manual or label. If in doubt, consult a specialist as this medical device may not be suitable for you.

## 6 attachments for whole body rehabilitation



You can work with all the major joints - wrist, elbow, shoulder, hip, knee and ankle, as well as trunk.



## Occupational therapy set and pelvic floor training adds-on



Interactive orthopaedic and EMG games make the rehabilitation process more interesting and fun for the patients. Automatic report availability.



*Motivation, enjoyment, and engagement are higher with the use of gamification [1]*

Luna EMG provides support for successfully facilitating older adults' physical activity through gamified technology.

(1) Kappen D.L., Mirza-Babaei P., Nacke L.E. (2020) Older Adults' Motivation for Physical Activity Using Gamified Technology: An Eight-Week Experimental Study. In: Gao Q., Zhou J. (eds) Human Aspects of IT for the Aged Population. Healthy and Active Aging. HCII 2020. Lecture Notes in Computer Science, vol 12208. Springer, Cham. [https://doi.org/10.1007/978-3-030-50249-2\\_22](https://doi.org/10.1007/978-3-030-50249-2_22)

## Training in all stages of rehabilitation

### MMT Levels

0 → 1 → 2 → 3 → 4 → 5



MMT LEVELS	EXPLANATION
0	No contraction
1	Flickering contraction
2	Full Range of Motion with eliminated gravity
3	Full Range of Motion with Against gravity
4	Full Range of Motion with Against gravity with minimal resistance
5	Full Range of Motion with Against gravity with maximal resistance

Getting to function - go step by step

From passive to active-resistive training

**LUNA EMG IS DESIGNED FOR A PATIENT'S COMPLETE RECOVERY JOURNEY**



*„Active participation of the patient contributes to significantly higher activation of the sensorimotor network during active motor control rather than during movement performed passively [2]”*

*„Robotics as an ideal means of training for severely affected patients where external assistance such as actuator assistance to movement and/or exoskeleton support may overcome problems of muscle weakness. [3]”*

## Research results




Article

### The Influence of EMG-Triggered Robotic Movement on Walking, Muscle Force and Spasticity after an Ischemic Stroke

Patrycja Lewandowska-Sroka <sup>1</sup>, Rafal Stabrawa <sup>1</sup>, Dominika Kozak <sup>2\*</sup>, Anna Pościwita <sup>2</sup>, Barbara Lyson-Ukłańska <sup>3</sup>, Katarzyna Bienias <sup>3</sup>, Anna Rokseła <sup>4</sup>, Marcin Kliś <sup>4</sup> and Michał Mikulski <sup>2</sup>

- Group:** 30 patients
- Training period:** 6 weeks, 1,5 h per day, 5 days a week
- Intervention:** Rehabilitation was based on individual standard physiotherapy and lower limb training with robot Luna EMG
- Results:** gait function improvement after the therapy



*“Due to a stroke I had a “frozen shoulder” which was very painful. While exercising actively on Luna EMG with shoulder extension the pain decreased and my range of motion in the shoulder increased. Now I can drive my car again and be more independent!”* Jacek, stroke patient

(2) Ziejka, K.; Skrzypek-Czerko, M.; Karłowicz, A. The Importance of Stroke Rehabilitation to Improve the Functional Status of Patients with Ischemic Stroke. *J. Neurol. Neurosurg. Nurs.* 2015, 4, 178–183, doi:10.15225/PNN.2015.4.4.6.  
 (3) Hesse, S. (2007). What does the clinician expect from machines and robots in Neurorehabilitation. In *European Symposium Technical Aids for Rehabilitation-TAR 2007*. Berlin (Germany).

