



REPAIR



REGENERATE



REPLENISH

# PRODUCT- PORTFOLIO



# CONTENTS

---

- 01 | INTRODUCTION
- 02 | WHAT IS MICRONEEDLING?
- 03 | INGREDIENTS
- 04 | INDICATIONS
- 05 | OUR PRODUCTS
- 06 | CONTACTS

# WE ARE THE INVENTORS OF THE MICRONEEDLING AND PIONEERS OF MODERN-DAY SKIN CONCEPTS

## Our Company Skin Care Mission & Vision

As the industry's trusted partner in microneedling, Dermaroller is committed to enhancing your practice.

With over 20 years of experience, we deliver cutting-edge, scientifically backed skin care solutions that prioritize patient safety and satisfaction. Join us on a transformative journey towards achieving remarkable results and empowering your clients' natural beauty.

## Revolutionizing The Market Since 1999

For 20+ years, Dermaroller has pioneered medical aesthetics with cutting-edge microneedling technology and advanced skin care solutions.

Trusted by professionals worldwide, our innovative approach sets new industry standards, delivering exceptional patient outcomes and redefining skin care practices.

## Research & Development

At Dermaroller, we are driven by a relentless pursuit of scientific innovation in the field of medical aesthetics. With a focus on research, development, and clinical expertise, we push the boundaries of skin care to deliver cutting-edge solutions.

Our commitment to excellence and transformative results sets us apart, making us the preferred choice for medical professionals worldwide. Discover the intersection of science and beauty as we revolutionize the future of skin care.

## Our Core Values

**Advancement:** Spearheading groundbreaking skin care solutions for medical professionals.

**Expertise:** Leading the way with cutting-edge skin care solutions for medical professionals.

**Efficacy:** Revolutionizing beauty with clinically-proven and long-lasting outcomes.

## PARTNERSHIP WITH DERMAROLLER GMBH



### Our Difference

Part of your choice is finding the right partner that best fits with your organisation and has the level of understanding and experience required to provide the impact your staff and clients will demand.



### Medical Excellence

We aim to be the partner of choice for our customers who themselves strive to deliver excellence in the medical field. We believe our technology & product's are the best in our field.



### Team

- ✔ Ensures we meet all regulations worldwide & safety is our no. 1 priority.
- ✔ Supports you in orders and after sales queries, marketing and customer support.
- ✔ Built and manufactured in Germany at our HQ with constant high quality checks.
- ✔ Customize all training requirements to the needs of your company.
- ✔ Ensures stock is available and delivered on time to you.



### Understanding

We have the pedigree, expertise and approach to make this project a success and provide a return on investment for you. Our understanding in your core strategy and development capability, makes Dermaroller the perfect choice of partner.



### Training

We aim to make a positive difference to practitioners by providing tools that will help them learn effectively and provide outstanding service to your clients.

### What is microneedling?

Microneedling also known as “Percutaneous Collagen Induction” (PCI) or “Collagen Induction Therapy” (CIT), is a skin treatment in which the outer layers of the skin, epidermis and dermis, are punctured in a minimally invasive manner, activating the body's own repair process and the renewal of damaged skin tissue by increasing collagen production.

With the CIT, the wound healing mechanisms are activated without causing side effects or long lasting damage to the epidermis; the puncture channels close again quickly.

### How does microneedling correct skin conditions?

Microneedling can be used to help patients with scarring, acne scarring, striae, burn scars and wrinkles / aged skin and, in conjunction with our XCellarisPRO skin care products, the entire treatment can reduce dry skin and hyper-pigmentation as well as the chances of hair thinning/loss.

### How does collagen induction therapy work?

When you perform microneedling on the skin at a depth as little as 0.5mm this releases your thrombocytes & erythrocytes which causes a spontaneous reaction, in turn releases your growth signals & factors.

The induced Angiogenesis is the physiological process through which new blood vessels form out of pre-existing vessels. This increases the blood flow back into the area of treatment and can promote growth.

The small pricking channels created react, and then we see an invasion of erythrocytes & leucocytes, due to improved cell to cell communication. Fibroblasts join the reaction, and we then start to see synthesis of Elastin and Collagen (Neocollagenesis).

This new matrix bolsters up the dermis from the inside and up towards the surface of the skin and closes the pricking channels. This is what we call tissue modulation.

The patient can continue this cascading effect from a professional microneedling treatment & at home by using the HC902 home care roller.



## Hero Ingredients in our products

### **MATRIXYL™ 3000**

is an anti-aging agent, approved in studies, which promotes production of collagen, and stimulates fibroblasts in the deeper layers of the skin. This then activates the growth of cells, in turn, improves texture and tone.

### **COPPER-TRIPEPTIDE-1**

is an anti-aging active for skin and hair, approved in studies! It regulates multiple regenerative and protective processes in the human body, including the skin.

### **HYALURONIC ACID**

improves the protective barrier and supports healing. It is an antioxidant: free radicals can break down collagen and elastin – antioxidants deactivate free radicals.

### **VITAMIN C**

is an antioxidant that deactivates free radicals and can reduce hyper pigmentation.

### **VITAMIN E**

is an antioxidant that deactivates free radicals and prevents damages caused by UV-rays, it also prevents inflammation and raises the tissue's ability to keep moisture.

### **VITAMIN A**

addresses any abnormalities in the fibroblast and increases the number of fibroblasts available to produce collagen & elastin.

## XCellarisPRO Skin care

Using our XCellarisPRO skin product range with microneedling provides the perfect aftercare for patients.

This will promote and support CIT and the skin's wound healing mechanism, resulting in longer lasting and improved results.

In addition to the ingredients listed we have a wide range of essential oils and extracts that are used in our products.

## We Care

We are Dermaroller GmbH. We invented microneedling. A pioneering spirit that drives and motivates our internationally operating, medium-sized family business to this day.

Our roots lie in the field of aesthetic medical technology. In addition, we specialize in dermatological, surgical, and aesthetic markets.

Our mission is to care for and protect the skin in all its wonderful facets.



## MEDICAL & AESTHETIC INDICATIONS



### Facial / Acne scarring

Our scar treatment significantly reduces the appearance of facial scars and acne scars. You can treat various types of scars from macular scars, rolling scar, box scar & ice pick scars. Feel confident again in your skin.

### Scarring

Our scar treatment significantly reduces the appearance of body scars and improves overall skin texture. Unlock a smoother, more even complexion and regain your confidence with our specialized solution.



### Striae Distensae

Say goodbye to stretch marks with our specialized treatment. By targeting and reducing the visibility of stretch marks, we help you achieve smoother and more even-toned skin, boosting your confidence.



## INDICATIONS IN COMBINATION WITH THE XCELLARISPRO LINE



### Anti-Aging

Our anti-aging treatment effectively reduces visible wrinkles and fine lines, for smoother, more youthful skin. Experience the rejuvenating effect and bring back the natural glow of your skin.

### Hyperpigmentation

Achieve a more even and radiant complexion with our hyperpigmentation treatment. By addressing uneven skin tone and reducing the appearance of dark spots, we reveal your skin's natural beauty and glow.



### Dry Skin / Rosacea

Revitalize and hydrate your dry skin with our nourishing treatment. Experience the intense moisture infusion that replenishes and restores your skin's natural balance, leaving it soft, supple, and revitalized.

### Hair Loss

Combat hair loss and stimulate hair growth with our effective treatment. By targeting hair follicles and promoting healthy hair growth, we help you regain thicker, fuller, and more voluminous hair.



\*All products to be applied to the skin post needling once the skin looks healthy and intact



## XCellarisPRO Twist microneedling device

Equipped with the **most powerful** handpiece of the market which allows a **minimal treatment time** and **optimal results**.

High stroke frequency: **150 insertions per second**, enables a quick treatment with high precision.

Powered device built for **consistency**.

**Highest comfort** during the treatment for the user and patient.

**Precise needle penetration** depth adjustment on the handpiece.

A single device can be used to treat a wide range of **medical and cosmetic indications**.

The needle count in the cartridge ensures **atraumatic penetration** and **minimizes the pain sensation**.

**Medical grade needles** manufactured to perfection through several consecutive manufacturing stages to ensure 360 degree symmetry.



### Technical Data

**Rated voltage:** 100 V - 240 V~

**Frequency of power supply:**  
50/60 Hz+

**Power consumption:**  
0.07 A (100V) - 0.03 A (240V)

**Working frequency:**  
Selectable from 50 to 150  
insertions per second

**Dimensions of control unit:**  
114 x 86 x 183 (mm)

**Weight of handpiece incl. cable:**  
ca. 260 g

**Control unit weight:** 825 g

XCELLARISPRO TWIST  
EXCELLENCE IN MICRONEEDLING –  
CRAFTED IN GERMANY



WHERE PRECISION MEETS INNOVATION

## CONSUMABLES FOR USING THE XCELLARISPRO TWIST



### CIT Hyal

Ensures the needle cartridge glides over the skin when performing professional microneedling.

**Ingredients:** Hyaluronic acid, Glycerin



### XCellarisPRO Twist Cartridges

Single use, individually packed and sterile cartridges. Includes a flowback safety system to avoid cross contamination. The needle length can be variably adjusted on the handpiece up to 2.5 mm. The needles are high quality, stainless steel microneedles, made in Germany.

## PRODUCTS



### Home Care Dermaroller (HC902)

With correct care usable for around 50 treatments.  
Dermaroller with 9 needle rows.  
162 microneedles with a length of 0.2 mm.

**Recommended use at home:** 2 – 3 times a week for optimal home care support for professional skin needling treatment



### Deep Cleanser

Daily skin cleanser that effectively removes skin impurities, dead cells and makeup, purifies clogged pores and optimally prepares the skin to absorb the ingredients of the XCellarisPRO products.

**Ingredients:** Copper-Tripeptide-1, Matrixyl™ 3000, Coco-Glucoside, Aloe Vera, Rice extract, Orange peel oil, Passion flower extract



### Moisturizer

Rich moisturizing cream of natural ingredients that retains moisture in the skin and strengthens the connective tissues to reduce lines and wrinkles.

**Ingredients:** Copper-Tripeptide-1, Matrixyl™ 3000, Vitamin A, Vitamin C, Vitamin E



### CC Cream

Concealing day cream with SPF 30 which allows you to control the desired skin tone. Mattifying, smoothes skin and improves skin's elasticity.

**Ingredients:** Copper-Tripeptide-1, Matrixyl™ 3000, Snail secret, UVA, UVB



### HA Treat

Deep hydration serums for the face, neck and décolleté. Restores hydration level of mature & dry skin, reduces inflammation and irritations, and prevents premature aging.

**Ingredients:** Copper-Tripeptide-1, Matrixyl™ 3000, Hyaluronic acid, Urea, Lactic Acid, Coco-Glucoside, Amino acids



### Cell Booster

Light moisturizing cream with collagen boosting ingredients, supports moisture retention, fights free radicals, supports formation of collagen, encourages a younger, wrinkle-free skin.

**Ingredients:** Copper-Tripeptide-1, Matrixyl™ 3000, Hyaluronic acid, Vitamin C, Argan Oil/Jojoba, Vitamin E, Niacinamide, Amino acids



## Cell Activator

Stimulating night cream, actively reduces wrinkles and improves elasticity of mature skin.

**Ingredients:** Copper-Tripeptide-1, Matrixyl™ 3000, Vitamin A, Rutin



## Cell Repair

Regenerating night serum with Retinol. Helps the skin to actively recover overnight, accelerates the regeneration processes, regulates the formation of calluses on the skin & can lighten age spots.

**Ingredients:** Copper-Tripeptide-1, Matrixyl™ 3000, Hyaluronic acid, Vitamin A, Vitamin E, Collagen & Elastin



## Hydration Mask

Intensive moisturizing sheet mask, gives a wonderful glow, soothing ingredients, promotes the maintenance and improvement of skin's own firmness, as well as rich and optimal hydration.

**Ingredients:** Copper-Tripeptide-1, Matrixyl™ 3000, Hyaluronic acid, Snow mushroom



## Bright Skin

Niacinamide (Vitamin B3) activates the cell renewal of the skin. Additionally its brightening effect supports a soft, tight and elastic skin.

**Ingredients:** Niacinamide (Vitamin B3), Vitamin A, Vitamin C, Vitamin E



## High Light Serum

Regulates hyperpigmentation in combination with the Bright Skin lotion. With the anti-oxidant effect wrinkles can be reduced as well as the skin elasticity increased, due to the enhanced collagen production.

**Ingredients:** Niacinamide (Vitamin B3) Plant extracts, Glycyrrhiza Glabra (root extract), Phyllanthus Emblica (fruit extract), Morus Alba (root extract)





## Pure HA Ampoules

Pure HA is a moisturising serum for the face, neck and décolleté, packaged in convenient, travel-friendly ampoules. 7 per pack.

**Ingredients:** Copper-Tripeptide-1, Matrixyl™ 3000, Hyaluronic acid, Urea

## Hair Lotion & Hair Treatment Set



Consists of one Hair Lotion, four Hair Rejuvenation Ampoules and four Hair Repair Ampoules.

Ensures the hair is firmly anchored in the root and stimulates the hair follicle, which can reduce hair loss

This creates the perfect base for the growth of strong and thick hair. Additionally, it strengthens hair roots, enhancing the overall health and resilience of the hair. The lotion is known to increase the thickness and density of hair, contributing to a fuller and more voluminous appearance.



It could be used on clients suffering from male and female androgenetic alopecia as well as alopecia areata, providing a comprehensive approach to managing different hair loss conditions.

**Ingredients:** Zinc sulfate, Azelaic acid, Copper peptides, Serenoa Serrulata fruit extract, Vitamin B6, Caffeine, Biotin, Niacinamide

## PRODUCTS FOR PROFESSIONALS ONLY



### Medical Dermaroller (MC4XX)

With correct care, usable for one treatment only. Dermaroller with 4 needle rows. 72 microneedles with a choice of needle lengths.


**Lengths:** 0.5 mm, 1.0 mm, 1.5 mm, 2.5 mm



### Medical Dermaroller (MC9XX)

With correct care, usable for one treatment only. Dermaroller with 9 needle rows. 162 microneedles with a choice of needle lengths.


**Lengths:** 0.5 mm, 1.0 mm, 1.5 mm, 2.5 mm



### Medical Dermaroller - Double (MC1825)

With correct care, usable for one treatment only. Dermaroller with 18 needle rows and 324 microneedles. Good for large surface area's of the body.

**Lengths:** 2.5 mm



- The original product Dermaroller - **Made in Germany**
- **Highest needle quality** made of stainless steel
- Various needle lengths for **optimal results**
- **Medical device Class IIa**
- **Individually sterilized** (gamma radiation sterilization)

mi.to.pharm GmbH  
Wilhelm-Mast-Str. 9  
38304 Wolfenbüttel  
Tel.: +49 (0)5331 710 8480  
Fax: +49 (0)5331 710 8489  
Mail: [info@dermaroller.com](mailto:info@dermaroller.com)  
Web: [www.dermaroller.com](http://www.dermaroller.com)

Inventors of microneedling  
and pioneers of  
modern-day skin concepts

